⚠ FIELD SAFETY NOTICE ⚠

Hyphen BioMed

BIOPHEN ARIXTRA[®] CONTROL PLASMA Ref 224001 – LOT #F1600528

Typographical error on the flyer:

Erroneous acceptance range for the C1 control

Dear customer,

Our traceability indicates that you have received this following product from Hyphen BioMed:

Product name	Product reference	Lot n°
BIOPHEN ARIXTRA [®] CONTROL PLASMA	224001	F1600528

We inform you that we have identified on the flyer, an error on the acceptance range of the C1 control:

Erroneous acceptance range of C1	<u>Correct</u> acceptance range of C1
0.3 – 0.56 µg/ml	<u>0.36</u> – 0.56 µg/ml

Potential consequence of the error

If the calibration curve is validated with a C1 control between 0.3 and 0.35 µg/ml, dosage values in the lower zone may be underestimated to 0.05 µg/ml.

Impact analysis

The maximum underestimation of 0.05µg/ml is not significant.

The potential risk does not affect the patient clinical setting.

Corrective actions

1. By the distributor (NODIA)

- Communicate the corrected flyer attached of 31.08.2017 to end users.
- Distribute the kits in stock with the new corrected flyer

2. By the end user

- Use the corrected flyer here enclosed

Cause and actions implemented

The internal investigation has identified a typographical error.

All measures are taken in Hyphen BioMed's Quality system to prevent any recurrence.

The French Competent Authority ANSM has been informed about this communication.

We would like to apologize for the inconvenience caused by this situation and remain at your disposal. For any question or information, do not hesitate to contact NODIA.

In order to meet regulatory requirements, can you please confirm receipt of this mail and provide

the enclosed corrected flyer to all concerned people in your laboratory?

Many thanks on beforehand.

Kind Regards,

6JG 01&+# - 6GCO

NODIA

Belgium/Luxembourg:

Scherpenhoek 10
B-2850 Boom
Belgium
Phone: +32 3 830 4644

The Netherlands:

Keurenplein 5
NL-1069 CB Amsterdam
The Netherlands
Phone: +31 20 610 80 06

www.nodia.com
info@nodia.com

HYPHEN BioMed

ZAC Neuville Université - 155 rue d'Eragry
95000 Neuville sur Oise (France)
Tel: +33.1.34.40.65.10 - Fax: +33.1.34.48.72.36
Web site: www.hyphen-biomed.com

ENGLISH / FRANÇAIS

BIOPHEN® ARIXTRA® CONTROL

Control plasmas for the quality control of Arixtra® measurements with anti-Xa method /

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