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READING TABLE

| TESTS | ACTIVE INGREDIENTS | QTY (mg/cup.) | REACTIONS/ENZYMES | RESULTS | |
|------------|------------------------------------|---------------|---|-----------------------------|------------------------------|
| | | | | NEGATIVE | POSITIVE |
| ONPG | 2-nitrophenyl-βD-galactopyranoside | 0.076 | βD-galactosidase (Ortho-NitroPhenyl-βD-Galactopyranosidase) | colorless | pale yellow to bright yellow |
| <u>LDC</u> | L-lysine | 0.46 | Lysine DeCarboxylase | yellow-green to bluish-grey | blue to blue-violet |
| <u>ODC</u> | L-ornithine | 0.28 | Ornithine DeCarboxylase | yellow-green to bluish-grey | blue to blue-violet |
| <u>URE</u> | urea | 0.32 | UREase | yellow | pink / pink-violet |
| CIT | trisodium citrate | 0.132 | CITrate utilization | yellow to yellow-green | green to blue |

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READING TABLE

| TESTS | ACTIVE INGREDIENTS | QTY (mg/cupule) | REACTIONS/ENZYMES | RESULTS | |
|------------|-----------------------------------|-----------------|---|-----------------------------|------------------------------|
| | | | | NEGATIVE | POSITIVE |
| ONPG | 2-Nitrophenyl-D-galactopyranoside | 0.076 | βD-Galactosidase (ortho-nitrophenyl-βD-galactopyranosidase) | Colorless | Pale yellow to bright yellow |
| <u>LDC</u> | L-Lysine | 0.46 | Lysine decarboxylase | Yellow-green to bluish-grey | Blue to blue-violet |
| <u>ODC</u> | L-Ornithine | 0.28 | Ornithine decarboxylase | Yellow-green to bluish-grey | Blue to blue-violet |
| <u>URE</u> | Urea | 0.32 | Urease | Yellow | Pink/Pink-violet |
| CIT | Trisodium citrate | 0.132 | Citrate utilization | Yellow | Green to blue |

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READING TABLE

| TESTS | ACTIVE INGREDIENTS | QTY (mg/cup.) | REACTIONS/ENZYMES | RESULTS | |
|-------|--------------------|---------------|---------------------|------------------------|---------------|
| | | | | NEGATIVE | POSITIVE |
| CIT | trisodium citrate | 0.132 | CITrate utilization | yellow to yellow-green | green to blue |

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Customer Relationship Center
bioMérieux Benelux

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