



User Name  
User Address Line 1  
User Address Line 2  
User Address Line 3

### **URGENT FIELD SAFETY NOTICE**

Product Name: ADEPT®  
Reference: FSN-21-005  
Action Type: Field Safety Notice on Medical Device – Read updated IFU  
Date: 18 October 2021  
Device Type: Metal on Metal Hip Resurfacing

Reason for this Field Safety Notice:

This FSN is to inform you that the indications and contraindications in the current instructions for use for these devices have been updated to reflect the post market surveillance (PMS) data presented at the last product approval which was completed in May 2021.

Root Cause:

The previous indications and contraindications were too general and were not specific to each product family. Whilst they were acceptable when they were originally used, clinical practise has evolved and the PMS data now available can be used to make them more specific and reflective of current clinical practice.

Revised indications for use:

The ADEPT® Hip Resurfacing System is intended for use in resurfacing hip arthroplasty for reduction or relief of pain and/or improved hip function in skeletally mature patients having Osteoarthritis.

The ADEPT® Hip Resurfacing System is intended for patients who, due to their relatively younger age or increased activity level, may not be suitable for traditional total hip replacements. Patient selection factors which should be considered include:

1. Ability and willingness of the patient to follow instructions.
2. Ability to commit to rehabilitation programme and attend regular follow-up as advised by health guidelines.
3. A good nutritional state of the patient and full skeletal maturity.
4. Patient bone density, osteopenia or family history of osteopenia including avascular necrosis (AVN) or osteonecrosis.
5. Patient age including need to regain mobility.

**Forever Active**



6. Patient femoral head size and specifically the slight increased risk of revision for the 46 mm diameter head.
7. Female patients who may have an increased risk of revision compared to males.

Revised contraindications for use:

1. Patients with active or suspected infection or sepsis in or around the hip joint.
2. Distant foci of infections (which may cause hematogenous spread to the implant site).
3. Osteoporosis, metabolic disorder, or a condition of marked bone loss which could impair bone formation.
4. Cases where there is inadequate neuromuscular status, vascular insufficiency, muscular atrophy, incompetent or deficient soft tissue surrounding the joint, neuromuscular disease, poor skin coverage around hip joint or prior implantation that cannot provide adequate support or fixation for the prosthesis which would compromise implant stability or post-operative recovery, making the procedure unjustifiable.
5. Skeletally immature patients.
6. Severe dysplasia.
7. Significant acetabular or femoral deformity that may preclude safe placement of components.
8. Bone stock which in the surgeons opinion is inadequate to support the device or where absorption is apparent on roentgenogram, this includes patients with severe osteopenia and a family history of severe osteoporosis or severe osteopenia, osteonecrosis or avascular necrosis (AVN) with >30% involvement of the femoral head (regardless of FICAT grade).

Note – In cases of questionable bone density, a Dual Energy X-Ray Absorptiometry (DEXA) scan may be necessary to assess inadequate bone stock.

9. Patients with multiple cysts of the femoral head (>1cm).
10. Patients with known moderate to severe renal insufficiency.
11. Patients who are immunosuppressed with diseases such as AIDS or persons receiving high doses of corticosteroids.
12. Patients with known or suspected metal sensitivity (e.g., jewellery).
13. Females of child bearing age or planning on becoming pregnant in the near term due to unknown effects on the foetus of metal ion release.

**Forever Active**



14. Metal sensitivity.

Additional information:

Additional information is provided in the IFU including:

The probability of revision at 10 years is:

<b>ADEPT® size range</b>	<b>Group</b>	<b>At risk at time point</b>	<b>Mean % cumulative revision rate at 10 year</b>	<b>95% confidence intervals</b>	
46mm to 58mm	Whole dataset (46 to 58mm heads, male and female)	584	4.5	3.5	5.8
	Males	518	4.0	3.0	5.3
	OA	545	4.4	3.4	5.8
<b>Higher risk groups:</b>					
46mm to 58mm	Females	64	8.6	4.5	15.5
46mm head only		63	9.2	4.7	17.0

Please read to current instructions for use which can be found at the following:

<https://info.matortho.com/reg/400-482.pdf>

The legal manufacturer for these devices is:

MatOrtho Limited  
 19/20 Mole Business Park  
 Randalls Road  
 Leatherhead  
 England  
 KT22 7BA  
 SRN: GB- MF- 000007614

The EU Authorised Representative for these devices is:

MatOrtho Ireland  
 The Black Church  
 Saint Mary's Place  
 Dublin  
 Ireland  
 D07 P4AX  
 SRN: IE-AR- 000006856

[AR@matortho.eu](mailto:AR@matortho.eu)

**Forever Active**



The EU Importer for these devices is:

MatOrtho Europe BV  
Tour & Taxis  
Havenlaan 86C  
Box 204  
1000 Brussels  
Belgium  
SRN: BE-IM-000007870

Dissemination of this FSN:

This FSN has been sent to all current users of the device.

Please pass this Field Safety Notice to all those who need to be aware of it and, if relevant, pass the FSN on to any another organization where the device may have been transferred

The relevant National Competent Authorities have been advised of the FSN.

All lots of the following devices:

Catalogue Number	Description	GMDN Code	EMDN Code	UDI-DI	Basic-UDI-DI
172-525	ADEPT® Metal On Metal Cup 52mm with 46mm bore	64684	P09080302	05055455507883	050554555000015Z
172-541	ADEPT® Metal On Metal Cup 54mm with 46mm bore	64684	P09080302	05055455507890	050554555000015Z
172-545	ADEPT® Metal On Metal Cup 54mm with 48mm bore	64684	P09080302	05055455507906	050554555000015Z
172-561	ADEPT® Metal On Metal Cup 56mm with 48mm bore	64684	P09080302	05055455507913	050554555000015Z
172-565	ADEPT® Metal On Metal Cup 56mm with 50mm bore	64684	P09080302	05055455507920	050554555000015Z
172-581	ADEPT® Metal On Metal Cup 58mm with 50mm bore	64684	P09080302	05055455507937	050554555000015Z
172-585	ADEPT® Metal On Metal Cup 58mm with 52mm bore	64684	P09080302	05055455507944	050554555000015Z

Forever Active

Catalogue Number	Description	GMDN Code	EMDN Code	UDI-DI	Basic-UDI-DI
172-601	ADEPT® Metal On Metal Cup 60mm with 52mm bore	64684	P09080302	05055455507951	050554555000015Z
172-605	ADEPT® Metal On Metal Cup 60mm with 54mm bore	64684	P09080302	05055455507968	050554555000015Z
172-621	ADEPT® Metal On Metal Cup 62mm with 54mm bore	64684	P09080302	05055455507975	050554555000015Z
172-625	ADEPT® Metal On Metal Cup 62mm with 56mm bore	64684	P09080302	05055455507982	050554555000015Z
172-641	ADEPT® Metal On Metal Cup 64mm with 56mm bore	64684	P09080302	05055455507999	050554555000015Z
172-645	ADEPT® Metal On Metal Cup 64mm with 58mm bore	64684	P09080302	05055455508002	050554555000015Z
172-661	ADEPT® Metal On Metal Cup 66mm with 58mm bore	64684	P09080302	05055455508019	050554555000015Z
181-146	ADEPT® Metal on Metal Resurfacing Head 46mm	64683	P09080402	05055455508064	0505545550000263
181-148	ADEPT® Metal on Metal Resurfacing Head 48mm	64683	P09080402	05055455508071	0505545550000263
181-150	ADEPT® Metal on Metal Resurfacing Head 50mm	64683	P09080402	05055455508088	0505545550000263
181-152	ADEPT® Metal on Metal Resurfacing Head 52mm	64683	P09080402	05055455508095	0505545550000263
181-154	ADEPT® Metal on Metal Resurfacing Head 54mm	64683	P09080402	05055455508101	0505545550000263
181-156	ADEPT® Metal on Metal Resurfacing Head 56mm	64683	P09080402	05055455508118	0505545550000263
181-158	ADEPT® Metal on Metal Resurfacing Head 58mm	64683	P09080402	05055455508125	0505545550000263



If you require further information, please contact your usual representative in the first place or customer services on [customer.services@matortho.com](mailto:customer.services@matortho.com) +44 (0)1372 366300.

...

Forever **Active**

MatOrtho Limited | 19/20 Mole Business Park | Randalls Road | Leatherhead | Surrey | KT22 7BA | United Kingdom

T: +44 (0)1372 224 200 | [info@MatOrtho.com](mailto:info@MatOrtho.com) | [www.MatOrtho.com](http://www.MatOrtho.com)

Company registered in England: 7323441



## Field Safety Notice Distributor/Importer Reply Form

<b>1. Field Safety Notice (FSN) information</b>	
FSN Reference number*	FSN-21-005
FSN Date*	Pre-filled by manufacturer
Product/ Device name*	ADEPT® Metal on Metal Hip Resurfacing
Product Code(s)	See Attached FSN

<b>2. Distributor/Importer Details</b>	
Company Name*	
Account Number	
Address*	
Shipping address if different to above	
Contact Name*	
Title or Function	
Telephone number*	
Email*	

<b>3. Return acknowledgement to Sender</b>	
Email	<a href="mailto:customer.services@matortho.com">customer.services@matortho.com</a>
Distributor/Importer Helpline	+44 (0)1372 366300
Postal Address	19/20 Mole Business Park Randalls Road Leatherhead, Surrey, KT22 7BA
Deadline for returning the Distributor/Importer reply form*	Pre-filled by manufacturer/sender/requester

<b>4. Distributors/Importers (Tick all that apply)</b>		
<input type="checkbox"/>	*I confirm the receipt, the reading and understanding of the Field Safety Notice.	Distributor/Importer to complete or enter N/A
<input type="checkbox"/>	I have identified customers that received or may have received this device	

**Forever Active**



<b>4. Distributors/Importers (Tick all that apply)</b>		
<input type="checkbox"/>	I have attached customer list	
<input type="checkbox"/>	I have informed the identified customers of this FSN	Date of communication:
<input type="checkbox"/>	I have received confirmation of reply from all identified customers	
Print Name*		Distributor/Importer print name here
Signature*		Distributor/Importer sign Here
Date *		

Mandatory fields are marked with \*

It is important that your organisation takes the actions detailed in the FSN and confirms that you have received the FSN.

Your organisation's reply is the evidence we need to monitor the progress of the corrective actions.

Forever **Active**

MatOrtho Limited | 19/20 Mole Business Park | Randalls Road | Leatherhead | Surrey | KT22 7BA | United Kingdom

T: +44 (0)1372 224 200 | info@MatOrtho.com | www.MatOrtho.com

Company registered in England: 7323441