

2023-12-21

URGENT FIELD SAFETY NOTICE

Manufacturer SRN: DE-MF-000020091

FSCA Reference: 946521 Hemoconcentrator assembled after its shelf-life had expired

FSN Type: New

Affected Product: Refer to Annex I List of affected products

Unique Device Identifier(s) (UDI-DI): 4037691057361, 4037691076041, 4037691149035, 4037691298351, 4037691563527

Affected Batch No.: Refer to Annex I List of affected products

For Attention of: Users of the medical device listed above

Dear valued customer,

Maquet Cardiopulmonary GmbH (MCP) would like to inform you with this letter about a corrective action for the above-mentioned Hemoconcentrator due to assembly after its shelf-life had expired.

The primary function of a hemoconcentrator is the elimination of excess water, electrolytes, and/or metabolites from the vascular space of a patient during cardiac surgery. The secondary function is retention of formed cellular elements (red blood cells, platelets, etc.) and albumin in the vascular system. The tertiary function is the removal excess fluid to be drawn from the interstitial space of the patient thereby preventing and/or minimizing peripheral edema.

A hemoconcentrator permits the retention of corpuscular blood components and plasma proteins while allowing the removal of excess plasma water thereby concentrating vascular volume. Unbound, low molecular weight solutes are removed from the vascular system with excess water in the process of hemoconcentration.

Problem description

The component "Hemoconcentrator" was assembled into the finished product (Hemoconcentrator tubing set) after its shelf-life had expired

After internal testing, it was determined that only one component lot was affected. Therefore, this Field Action is limited to products that containing Hemoconcentrator of the one affected batch.

Hazardous situation

In course of a Health Hazard Evaluation (HHE), Maquet Cardiopulmonary GmbH determined the following hazardous situations for the expired Hemoconcentrators:

- Patient is exposed to inappropriate high Thrombogenicity
- Patient is exposed to inflammatory agents
- Patient is exposed to inappropriately high hemodilution

Potential harm

The possible immediate and/or long-range health consequences and risk levels of the non-conformance include the following:

- Coagulation disorder
- Ischemia (Thromboembolism)
- Bleeding
- Inflammation
- Anemia
- Hemodilution

Maquet Cardiopulmonary GmbH has not identified any complaints of patient harm, serious injuries, or deaths due to the failure modes described above.

FIELD SAFETY NOTICE



DMS No.: 3287630 V 01

Page: 4 of 9

We sincerely apologize for any inconvenience this may cause you and will do our utmost to carry through this action as swiftly as possible.

As required, we have provided this notification to the necessary Regulatory Agencies.

Should you have questions or require additional information, please contact your local Getinge representative, or send an e-mail to FSCA.cp@getinge.com.

Sincerely,

Contact details of manufacturer

Tom Peters
Maquet Cardiopulmonary GmbH
Kehler Str. 31
76437 Rastatt
GERMANY
Phone: +49 7222 932 - 0
Email: FSCA.cp@getinge.com

CUSTOMER RESPONSE FORM

FSCA Reference: 946521 Hemoconcentrator assembled after its shelf-life had expired

Affected Product: Refer to Annex I List of affected products

Affected Batch No.: Refer to Annex I List of affected products

Please send this form at the latest by **January 17, 2024**, to your local Getinge representative.

By completing this document and signing it, I acknowledge that I have read and understand the following associated points:

- I have read and understand this Field Safety Notice for affected product Hemoconcentrator. We will take action as soon as possible according to given instructions.
 - I confirm that I have distributed this Field Safety Notice to the affected personal.
- I do not have any Hemoconcentrator in my inventory.
- I have following Hemoconcentrator in my inventory.

| Article No. | Description | Batch No. | Quantity |
|-------------|-------------|-----------|----------|
| | | | |
| | | | |
| | | | |
| | | | |

Your Comments:

Country

Hospital / Clinic (full address)

Date

Name (Function)

Signature

FIELD SAFETY NOTICE



DMS No.: 3287630 V 01

Page: 6 of 9

Please return the completed form to your local Getinge representative by email enter local Getinge mail address or via post enter local Getinge address or FAX.

Template: CP-SOP-001-T-02 V02, Effective date 2019-09-15

Annex I List of affected products

This Annex I List of affected products is considered a supplementary attachment to the 946521 Field Safety Notice.

Canada:

| Article no. | Item Description | Batch no. |
|-------------|---|--------------------------|
| 701027710 | H 52570#BC 140 Plus with Lines | 3000273187 |
| 701035298 | BO-H 48771#BC 140 Hemoconcentrator Plus | 3000276250 3000277066 |

Czech Republic:

| Article no. | Item Description | Batch no. |
|-------------|--|------------|
| 701005142 | P-0400#Hemoconcentrator set incl. BC 140 | 3000273199 |

France:

| Article no. | Item Description | Batch no. |
|-------------|--|------------|
| 701005142 | P-0400#Hemoconcentrator set incl. BC 140 | 3000242892 |

Germany:

| Article no. | Item Description | Batch no. |
|-------------|--|--------------------------|
| 701005142 | P-0400#Hämokonzentrator Set mit BC 140 | 3000242892 3000273199 |

Netherlands:

| Article no. | Item Description | Batch no. |
|-------------|-------------------------------------|------------|
| 701048516 | H 31471#Hemoconcentrator Set BC 140 | 3000247617 |
| 701030315 | H 15981#Set BC 140 Plus | 3000247619 |

FIELD SAFETY NOTICE

DMS No.: 3287630 V 01

Page: 8 of 9

Philippines:

| Article no. | Item Description | Batch no. |
|-------------|---|------------|
| 701005142 | P-0400#Hemoconcentrator set incl. BC140 | 3000273199 |

Portugal:

| Article no. | Item Description | Batch no. |
|-------------|--|------------|
| 701005142 | P-0400#Hemofiltration Set incl BC140Plus | 3000242892 |

South Africa:

| Article no. | Item Description | Batch no. |
|-------------|---|------------|
| 701005142 | P-0400#Hemoconcentrator set incl. BC140 | 3000273199 |

Switzerland:

| Article no. | Item Description | Batch no. |
|-------------|--|--------------------------|
| 701005142 | P-0400#Hämokonzentrator Set mit BC 140 | 3000242892 3000273199 |

United Arab Emirates:

| Article no. | Item Description | Batch no. |
|-------------|--|------------|
| 701005142 | P-0400#Hemoconcentrator set incl. BC 140 | 3000242892 |

Annex II Further information regarding Hazardous situation, Harms and Risk Levels

This Annex II Further information regarding Hazardous situation, Harms and Risk Levels is considered a supplementary attachment to the 946521 Field Safety Notice.

| Hazardous situation | Harm | S from part III | P from above | Risk | | |
|--|---|-----------------------|--------------------|--------------------------|-------------------------------------|--------------------------|
| | | | | Low | Med | High |
| Patient is exposed to inappropriate high Thrombogenicity | Coagulation disorder ^{b,c} | 3 | 3 | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> |
| | Ischemia (Thromboembolism) ^b | 4 | 3 | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> |
| | Bleeding ^c | 3 | 3 | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> |
| Patient is exposed to inflammatory agents | Inflammation | 3 | 3 | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> |
| Patient is exposed to inappropriately high hemodilution | Anemia | 3 | 3 | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> |
| | Hemodilution | 3 | 3 | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> |
| Product exchange/replacement | User inconvenience | 2 | 2 | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> |

Severity Definitions:

Negligible (1) Inconvenience or temporary discomfort of patient, user or third party. No medical intervention or follow-up treatment is required

Low (2) Temporary injury or disability of patients, users or third parties. No medical intervention or follow up treatment is required.

Critical (3) Temporary injury or disability of patients, users or third parties. Medical intervention or follow-up treatment is required.

Catastrophic (4) Permanent injury or disability (e.g., loss of a body part), a life-threatening situation or death of patients, users or third parties

Probability Definitions:

Improbable (1) Harm is not likely.

Remote (2) Harm occurs infrequently

Occasional (3) Harm may occur occasionally / intermittent

Probable (4) Harm may occur often

Frequent (5) Harm will occur repeatedly