



Facility  
Service  
Address  
Address  
ZipCode City  
Country

# URGENT: FIELD SAFETY NOTICE

## Medical Devices Emergency Treatment

24/04/2026

ATTENTION: Pharmacist/Risk Manager responsible for medical device vigilance and the Biomedical/Engineering Department

### Field Safety Notice regarding the Emergency Treatment of Namic White Star Handle Manifolds included in Medline's Namic Systems

Medline Reference: FSN-26/06  
MoH Reference: N/A  
Product description: Namic Systems with Namic White Star Handle Manifolds included  
Legal Manufacturer SRN: FR-PR-000003713  
Action type: Field Safety Notice – Emergency Treatment  
Product codes: See Appendix 1

Dear Customer,

UPDATE 24<sup>th</sup> APRIL 2026

This letter follows up on the Field Safety Notice – Emergency Treatment regarding the Recall of Namic White Star Handle Manifolds issued on 21<sup>st</sup> April 2026. We need to inform you that additional product codes fall under the scope of the Field Safety Notice (FSN) of Namic White Star Handle Manifolds. All impacted product codes for this FSN are listed in Appendix 1 (page 5) attached hereto.

An additional required action has been added to this FSN to ensure continuity of use. Customers are now instructed to order a replacement single sterile manifold to replace the defective manifold included in the NAMIC Kit.



Medline hereby informs all customers in the Netherlands that single sterile replacement manifolds are available upon request to replace the defective Namic White Star Handle Manifolds included in NAMIC Systems.

To ensure continuity of interventional cardiology procedures, customers are instructed to order a replacement single sterile manifold to replace the defective manifold supplied within the NAMIC Kit.

In exceptional circumstances where an emergency procedure must be performed, no alternative device is available, and postponing the procedure would place the patient at significant risk, the use of the originally supplied Namic White Star Handle Manifold may be considered as a last resort.

In such cases, the user must perform a full and careful flush of the Namic White Star Handle Manifold and associated tubing with saline, thoroughly inspect the manifold and lines for any particulate matter, and must not use the device if any particulate is observed.

**REASON FOR THE FIELD SAFETY NOTICE:**

Customers in the Netherlands have been disproportionately impacted by the limited availability of Namic White Star Handle Manifolds following previous corrective actions.


This Field Safety Notice is issued to inform customers that single sterile replacement manifolds are available upon request to replace the defective manifolds supplied within NAMIC Systems, and to define controlled measures to ensure that patients at immediate and significant clinical risk may continue to receive emergency interventional cardiology treatment.

**POTENTIAL RISKS:**

The potential presence of foreign particulates within the fluid pathway of Namic White Star Handle Manifolds may, under certain conditions, result in:

- Delay in patient care
- Embolism
- In rare instances, death

Based on review of currently available post-market surveillance data, there have been no reported events resulting in embolism or death.





**ACTIONS REQUIRED:**

Step 1: Please take note of this FSN and inform all users in your facility.

Step 2: Please complete the Acknowledgment Receipt (*page 4*). Then, return page 4 by email as soon as possible **but no later than 15<sup>th</sup> May 2026.**



Step 3: To obtain a single sterile manifold as a replacement for the defective component identified in your NAMIC Kit, please contact your Medline Sales Representative to place an order.

Step 4: If you have any questions or need support, please contact your Medline Sales Representative or your local Medline Customer Service.

The relevant competent authorities have been informed of this FSN.

Please proceed with the following pages to acknowledge receipt of this FSN. Medline would like to thank you for your cooperation.

**WARNING LABEL:**

	<b>FSN-26/06</b>	
<b>DON'T USE the Namic White Star Handle Manifold included in Medline's Namic Systems, except in life-threatening emergencies.</b>		
<b>If so, please FLUSH and INSPECT for loose particles before use</b>		

Yours sincerely,

[Redacted signature]

*This urgent and important field safety notice is only addressed to Dutch facilities that may receive the referenced products.*

[Redacted contact information]

[Redacted contact information]



**Please email the Acknowledgement Receipt to the following email address:  
GMB-EU-FSN-FSCA-KLEVE@medline.com**

**Medline Reference: FSN-26/06**

Please complete the Acknowledgement Receipt, and send back by email as soon as possible, but no later than 15<sup>th</sup> of May 2026. The referenced products are listed in the Appendix (page 5).

By completing and signing this acknowledgement receipt, I confirm that I have read and I understood the instructions provided. I acknowledge receipt of the FSN-26/06 by signing this document and returning it to Medline.

This notice needs to be passed on to all those who need to be aware within your organization or to any organization where the potentially affected devices have been transferred. (As appropriate)

Please transfer this notice to other organizations on which this action has an impact. (As appropriate)

Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.

Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback.

If you are a dealer, wholesaler, distributor/reseller, that distributed any referenced products to other facilities: per Medical Device Regulation 2017/745, Article 14, part 4, please distribute this notification to your customers and provide confirmation to Medline that your customers have been notified by completing the information below and returning it to Medline per email to [REDACTED].

Date: \_\_\_\_\_

Name: \_\_\_\_\_

Position: \_\_\_\_\_

Facility or Business Entity: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

Medline Account Number: \_\_\_\_\_

Telephone: \_\_\_\_\_

Email address: \_\_\_\_\_

Signature: \_\_\_\_\_

**Medline International France SAS**



**Quality & Regulatory Affairs Dept.**





## Appendix 1

Item number
606000776
606209316A
60610625918
60148473821
606209331
6060007128B
606209318
NABK60008
NADK60002
NANK60013
NASK60002
NASK60003
NAUK60003
NAUK60006A
NAWK60001
NAWK60002

Item number
NAWK60004
NAWK60006
NAWK60008
NAWK60012
NAWK60016
NAWK60019
NAWK60020
NAWK60021
NAWK60022
NAWK60023
NAWK60025
NAWK60027
NAWK60028
NAWK60029
NAWK60036