Alle Anwender (Betreiber, Vertreiber, Patienten etc.)
von N-LINEpro mit der Software-Version: xx_2.0.39/2.01
All users (operators, distributors, patients, etc.)
by N-LINEpro with the software version: xx_2.0.39/2.01

Herford, den 08.07.2025

Dringende Sicherheitsinformation / Field Safety Notice

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In the interest of patient safety we would like to inform you about the following safety measures.

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The following devices in connection with the referenced input software are concerned:

<u>Artikel-Nr.</u>	<u>Modellbezeichnung</u>	<u>Article No.</u>	<u>Model Name</u>
4557	N-LINEpro UVB 311nm	4557	N-LINEpro UVB 311nm
4558	N-LINEpro UVA	4558	N-LINEpro UVA
4559	N-LINEpro Universale	4559	N-LINEpro Universal

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You can find the serial number on the product's nameplate. The nameplate is located in the middle of the inner, left-hand entrance area of the cabin. The software version is displayed in the info menu of the cabin.

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MEDlight GmbH has received information about an individual case in which a very fast and very unusual operating sequence, the entry of an mJ value instead of a J value, the non-observance of visual and acoustic warnings regarding both the dose and the duration of treatment (manual

override of a soft stop by the user) and the non-monitoring of the application by the user could have led to an undesired dose exposure of the patient, in which damage to the patient's health due to burns cannot be ruled out according to the current state of knowledge.

According to the examination results of MEDlight GmbH, this can only occur in very rare cases and only after the sequence of events described below in combination with the described non-observance of the visual and acoustic warnings by the user:

- The user enters a value in mJ/cm² instead of a value in J/cm² in the field for entering the dosage.
- Shortly after entering the value in mJ/cm², the user jumps back to the main menu and immediately (within a few milliseconds) back to the treatment parameterization, while the initial input in mJ/cm² is still being processed in the background.
- In the treatment parameterization, the user does not wait for the "Start" button to appear, but presses the space on the touchscreen where the "Start" button normally appears, even before or at the exact moment when the "Start" button is displayed.
- By entering a dosage value in mJ/cm², the very fast jumps between the treatment parameterization and the main menu as well as the very fast triggering of the start without waiting for the "Start" button to appear, the system may make the error that the last entered mJ value is interpreted as a value in the unit J/cm² and set as the treatment dose.
- Regardless of this, even in this rare case, the system correctly triggers both a visual and an acoustic warning ("Warning level exceeded") before the start of treatment, alerting the user to the high dose and prompting them to check and ensure that the selected dose is correct (soft stop).
- In addition, in this rare case, the system actively indicates a very long treatment time of several hours, which is not plausible (in the specific individual case that occurred: 9 hours).
- Only if the user then ignores both the visual and acoustic warning of the high dose and the indication of the implausibly long treatment time according to the operating sequence described above, which is not intended in this way, and still triggers the start of treatment manually, can treatment begin with a dose in J/cm² that was originally entered in mJ/cm².
- If then the user does not comply with the requirement to monitor the treatment set forth in the instructions for use, this can lead to an overdose of UV irradiation to the patient with corresponding damage to health, including burns.
- According to the findings of MEDlight GmbH, this very rare error scenario requires exactly the above-mentioned and precisely timed interaction between an mJ/cm² input, very fast mask changes and a very fast "Start" click in combination with the above-mentioned skipping of the various warnings on dose and treatment duration and, according to current knowledge, can only be reproduced under these special circumstances in individual cases. According to the current knowledge of MEDlight GmbH, the error has not yet occurred in other cases.

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Potential danger and possible risks for patients:

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In the event of the above-described unintended sequence of operations, and if acoustic and visual warnings are ignored, therapy may, in very rare cases, be initiated with an excessive dose by a negligent user. This can lead to the treatment of patients with excessive doses and, especially if the treatment is not monitored, to burns.

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Until new input software is available, the dose should only be entered in the software in J/cm^2 . The user should refrain from selecting doses in the unit range J/cm^2 as described in the user manual under 3.3.4 “Selecting dose units.” Warning notices must be observed and checked.

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MEDlight will provide a software update for the affected products, which will prevent treatment from starting even if the warning messages are ignored in situations similar to the one described above (hard stop). In addition, we are investigating whether it will be possible to only allow dose entries in J/cm^2 in future. No further measures are necessary.

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Please ensure in your organization that all users and potential users of the above products and other persons to be informed are aware of this “Field Safety Notice”.

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If you have redistributed the products which are subject to this field safety corrective action, please identify the affected customers and inform these customers immediately of this field safety corrective action by forwarding a copy of this “Field Safety Notice” or informing the contact person indicated below. We kindly ask you to monitor and follow-up on this measure at your customers!

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Please complete the enclosed reply form in full and send it to the specified address (by mail, e-mail or fax) on or before July 31, 2025. Please return the reply form even if you have no remaining stock. Due to legal requirements, MEDlight GmbH must ensure and document that you have received this information and your reply serves as evidence.

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The Federal Institute for Drugs and Medical Devices (BfArM) has received a copy of this “Field Safety Notice”.

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If you have any questions, please do not hesitate to contact us.

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Confidential business information

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We apologize for any inconvenience caused to you in connection with the implementation of this field safety action and would like to thank you in advance for your support in order to enable a quick and effective implementation of this measure in the interests of patient and user safety.

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Annex to “Field Safety Notice”

Reply form

Upon receipt of your response form, you will be contacted by our customer service to arrange a software update of the concerned products. Only the input software of the products with concerned serial numbers will be updated.

Please do not return any products unsolicitedly!

Please tick the appropriate fields!

☐ **Operator Organization**

☐ We have distributed this “Field Safety Notice” from July 8, 2025 to our users.

☐ We and our users have read and understood this “Field Safety Notice” from July 8, 2025.

☐ We have equipped the affected products with the safety notice sent to us.

(Please attach a list of the serial numbers/UDIs of the affected products!)

☐ **Distributor**

☐ We have read and understood this “Field Safety Notice” from July 8, 2025.

☐ We have informed our customers by forwarding this “Field Safety Notice” (Please attach a list of affected serial numbers/UDIs located at your customers)

☐ We do not have **any remaining stock** of the affected products.

☐ We have the **following quantities** in stock: _____
(Please attach a list of the serial numbers/UDIs of the affected products!)

Please enter your contact details here:

Name of the organization: _____

Postal address: _____

Name of the contact person: _____

Telephone & Fax: _____

E-Mail: _____

Date

Signature